

Stabilizer ROM & Fixed Ankle Walking Boot

ROM (Range of Motion) Joint Features :

- The range of motion joint can be fixed in any position between 45° of dorsi-flexion and 45° of plantar-flexion
- The range of motion joint adjusts in 7.5° increments in both dorsi-flexion and plantar-flexion, and can either be completely locked out at any desired interval or left with desired range.



ROM Walker

Stabilizer Features :

- Light weight and low profile construction
- Non-skid rocker bottom sole evenly progresses the patient from heel to toe
- Fixed malleable uprights provides custom fit ensuring stability
- Integrated strapping system provides full circumferential compression of the limb increasing stability and comfort
- Fixed malleable uprights provides custom fit ensuring stability
- Sturdy injection molded footpiece with padded heel and footplate
- Breathable and washable fiber liner with sewn-in heel pad for added comfort



Fixed Ankle Walker

Indications and Uses :

- Acute ankle sprains
- Soft tissue injuries of the lower leg
- Ankle stress fractures
- Stable fractures of the ankle and lower leg
- Stress and stable fractures of the foot
- Achilles tendonitis
- Achilles tendon repairs

ALSO AVAILABLE WITH INFLATABLE AIR LINER



Air Pump Liner

P/N (with air liner)	Description	P/N (no air liner)
09-1010	Range of Motion Walker, Small	09-1018
09-1012	Range of Motion Walker, Medium	09-1019
09-1013	Range of Motion Walker, Large	09-1020
09-1001	Fixed Ankle Walker, Small	09-1021
09-1002	Fixed Ankle Walker, Medium	09-1022
09-1003	Fixed Ankle Walker, Large	09-1023
09-1015	Replacement Liners, Small	09-2020
09-1016	Replacement Liners, Medium	09-2022
09-1017	Replacement Liners, Large	09-2023