

281

Ligaflex Boxer and Ligaflex Finger Braces



b inst s.

1 Fal



Ligaflex[®] Boxer

Wrist-Hand-Finger Orthosis PDAC Approved L3807/L3809

Indications

- Non-displaced 4th and 5th metacarpal fracture
- Postoperative 4th and 5th metacarpal support
- Following strict cast immobilization

Comfort

- Soft fabric with foam liner to cover the injured area
- 3D inner breathable fabric for daily comfortable wear
- Finger fabric separator to protect the immobilized finger

Immobilization

- Two shapable rigid stays for a strict immobilization from the wrist up to the finger to aid in healing
- Pre-shaped palmar and dorsal stays (dorsal = 90° and palmar = 20°) for optimum immobilization

User-Friendly

- Full opening of the brace for easier application and positioning
- Convenient buckle loop system to tighten the brace
- Machine washable

Ligaflex[®] Finger

Wrist-Finger Orthosis Recommended PDAC L3923/L3924

Indications

- Long finger injuries including sprains, dislocations and fractures (functional treatment or post surgical support)
- After-effects of injuries (sprain, dislocation and fracture)

Comfort

- Four straps, two fixed and two removable, for optimized support and fit to improve wearing comfort
- Wrist strap with anatomical shape to improve wearing comfort
- Curved top of the stay to cover and protect the finger

Immobilization

- Unique stay with four straps to stabilize the brace
- Rigidity for strict immobilization from the wrist up to the finger
- Adaptability to fit the patient's wrist and hand shape and follow the finger angulation
- Two stay positions, palmar and dorsal, for optimum immobilization

User-Friendly

- Full opening of the brace for easier application and positioning
- Machine washable







Ligaflex[®] Boxer Wrist-Hand-Finger Orthosis

Size	Wrist Circumference	
1	5 — 5 ³ /4"	~
2	6 — 6 ³ /4"	
3	7 — 7 ³ /4"	
4	8 - 9"	V -

Ligaflex[®] Finger

Wrist-Finger Orthosis

Size	Tip of Injured Finger	
1	< 5"	\wedge
2	5 - 6 1/4"	T/P
3	6 1/4 - 7 1/2"	
4	7 ¹ / ₂ - 8 ³ / ₄ "	

5545 St. Jacques O. Montreal, Qc H4A 2E3 O: 514.369.3311 Toll Free / Sans Frais : 1.877.222.3311 🍁 Fax Toll Free / Telec Sans Frais : 1.877.527.1911 www.bmiortho.com