

# **ETHNOCARE**

---

# **LINER**



**PRODUCT CATALOG**

# Liner size

Silicone liner			
Liner size	Length	Circumference 4cm above distal end	Circumference proximal end
16	35 ± 2 cm	16	20
18	35 ± 2 cm	18	21.5
20	35 ± 2 cm	20	23.5
21	35 ± 2 cm	21	24.5
22	35 ± 2 cm	22	25.5
23.5	35 ± 2 cm	23.5	27
25	35 ± 2 cm	25	28.5
26.5	35 ± 2 cm	26.5	30
28	35 ± 2 cm	28	31.5
30	35 ± 2 cm	30	33.5
32	35 ± 2 cm	32	35.5
34	35 ± 2 cm	34	37.5
36	35 ± 2 cm	36	39.5
38	35 ± 2 cm	38	41.5
40	35 ± 2 cm	40	43.5
42	35 ± 2 cm	42	45.5
45	35 ± 2 cm	45	48.5
50	35 ± 2 cm	50	58.5
55	35 ± 2 cm	55	63

Gel			
Liner size	Length	Circumference 4cm above distal end	Circumference proximal end
S	40 ± 2 cm	15-22	24-35
M	40 ± 2 cm	19-29	26-40
L	40 ± 2 cm	23-27	37-61
XL	40 ± 2 cm	27-40	40-64
XXL	40 ± 2 cm	31-53	42-70

# TT Locking Silicone Liner

Made from soft, medical-grade silicone, the TT Locking Silicone Liner provides lasting comfort and gentle support. Its seamless textile construction allows easy donning and removal, while the integrated reinforcement minimizes distal stretch for a stable fit. The distal end features an anti-rotation connector compatible with all standard pin systems.

## Features



Soft  
Silicone



Fabric  
cover



Integrated  
Matrix



Distal  
locking




## Specs

Product details		
Wall thickness	3 mm or 6 mm	
Matrix	Integrated in knitted fabric	
Distal connection	With	
Color:	Light grey	
Size	16, 18, 20, 21, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 38, 40, 42, 45	
Item No		
3 mm	LNR-SIL-L-3-XX	(xx = Size )
6 mm	LNR-SIL-L-6-XX	(xx = Size )

## Size Determination:

Measure 4 cm above the distal end of the residual limb and select the corresponding or next smaller size.



3mm profile / 6mm profile

**ETHNOCARE**

# TT Cushion Silicone Liner

The TT Cushion Silicone Liner combines softness and durability for all-day comfort. Made from medical-grade silicone, it gently conforms to the limb for even pressure distribution. Its seamless textile exterior ensures smooth rolling and easy donning, while the cushioning design provides secure suspension without the need for a pin locking system.

## Features



Soft  
Silicone



Fabric  
cover



Integrated  
Matrix




## Specs

Product details		
Wall thickness	3 mm or 6 mm	
Matrix	Integrated in knitted fabric	
Distal connection	Without	
Color:	Light grey	
Size	16, 18, 20, 21, 22, 23.5, 25, 26.5, 28, 30, 32, 34, 36, 38, 40, 42, 45	
Item No		
3 mm	LNR-SIL-C-3-XX	(xx = Size )
6 mm	LNR-SIL-C-6-XX	(xx = Size )

## Size Determination:

Measure 4 cm above the distal end of the residual limb and select the corresponding or next smaller size.



3mm profile / 6mm profile

# TF Locking Silicone Liner

Made from soft, medical-grade silicone, the TF Locking Silicone Liner provides lasting comfort and gentle support. Its seamless textile construction allows easy donning and removal, while the integrated reinforcement minimizes distal stretch for a stable fit. The distal end features an anti-rotation connector compatible with all standard pin systems.

## Features



Soft  
Silicone



Fabric  
cover



Integrated  
Matrix



Distal  
locking




## Specs

Product details	
Wall thickness	2 mm
Matrix	Integrated in knitted fabric
Distal connection	With
Color:	Light grey
Size	25, 26.5, 28, 30, 32, 34, 36, 38, 40, 42, 45, 50, 55
Item No	
2 mm	LNR-TF-SIL-L-2-XX (xx = Size )

## Size Determination:

Measure 4 cm above the distal end of the residual limb and select the corresponding or next smaller size.



# Gel Locking Liner

Made from soft, medical-grade gel, the Gel Locking Liner delivers exceptional comfort and cushioning. Its functional textile exterior ensures smooth rolling for easy donning and removal. The distal connector provides a secure attachment and is compatible with all standard pin systems.

## Features



Soft  
Gel



Fabric  
cover



Distal  
locking




## Specs

Product details					
Wall thickness	3 mm, 6 mm or 9mm (Back / Front)				
Matrix	Without				
Distal connection	With				
Color	Dark grey				
Size	S	M	L	XL	XXL
Circumference 4cm above distal end in cm	(15-22)	(19-29)	(23-37)	(27-40)	(31-53)
Circumference at proximal end in cm	(24-35)	(26-40)	(37-61)	(40-64)	(42-70)
Item No					
3/3 mm	LNR-GEL-L-3-XX		(xx = Size )		
3/6 mm	LNR-GEL-L-3-6-XX		(xx = Size )		
3/9 mm	LNR-GEL-L-3-9-XX		(xx = Size )		

## Size Determination:

Measure 4 cm above the distal end of the residual limb and select the corresponding or next smaller size.



**ETHNOCARE**

# Gel Cushion Liner

Crafted from soft, medical-grade gel, the Gel Cushion Liner provides lasting comfort and gentle cushioning throughout the day. Its smooth, functional textile makes donning and removal easy, while the pin-free design offers a secure, comfortable fit for users who prefer suction or sleeve suspension systems.

## Features



Soft  
Gel



Fabric  
cover




## Specs

Product details					
Wall thickness	3 mm, 6 mm or 9mm (Back / Front)				
Matrix	Without				
Distal connection	Without				
Color	Dark grey				
Size	S	M	L	XL	XXL
Circumference 4cm above distal end in cm	(15-22)	(19-29)	(23-37)	(27-40)	(31-53)
Circumference at proximal end in cm	(24-35)	(26-40)	(37-61)	(40-64)	(42-70)
Item No					
3/3 mm	LNR-GEL-C-3-XX		(xx = Size )		
3/6 mm	LNR-GEL-C-3-6-XX		(xx = Size )		
3/9 mm	LNR-GEL-C-3-9-XX		(xx = Size )		

## Size Determination:

Measure 4 cm above the distal end of the residual limb and select the corresponding or next smaller size.



**ETHNOCARE**

## Our products :



### OVERLAY TT

The Overlay TT is a lightweight, breathable sleeve that provides uniform pressure and support inside the socket. Its air expansion system adapts to different limb shapes for a customized fit. Easy to wear under clothing, it enhances comfort throughout the day.



### OVERLAY TF

The Overlay TF offers transfemoral users a comfortable, adjustable fit with consistent pressure distribution. Its advanced textile and air system adapt to the residual limb for improved suspension and comfort. Slim and unobtrusive, it works seamlessly with most sockets.



### UNDERLAY

The Underlay is a moisture-wicking sleeve worn under the prosthetic liner to keep the limb cool and dry. Its special knit design absorbs sweat and releases it out of the liner, reducing discomfort. It ensures lasting freshness throughout the day.

#### Contact us :



Ethnocare.ca




clinics@ethnocare.ca



+1 (418) 934-5669



240 rue saint Jacques-O, 501,  
Montréal, QC, Canada, H2Y 1L9



Website